



MRI GE DISCOVERY MR750 3T - CR CHU SAINTE-JUSTINE					
RESTING STATE FMRI ACQUISITION PROTOCOL					
T2-WEIGHTED					
SEQUENCE NAME	CUBE T2W	TE	160	ECHO TRAIN LENGHT	200
TYPE OF SEQUENCE	3D Isotropic Fast SE	TR	5000	FREQ. DIR	S/I
IMAGE ENHANCEMENT	PURE FILTER	FLIP ANGLE	90	PIXEL BW	195.312
				SLICE THICKNESS	0.75
RS-fMRI (AP)					
SEQUENCE NAME	GRE-EPI	TE	30	#VOLUMES	100
TYPE OF SEQUENCE	Single-shot GRE-EPI	TR	3000	FREQ. DIR	R/L
SLICE ORDER	Interleaved bottom/up	FLIP ANGLE	70	PIXEL BW	7812.5
				SLICE THICKNESS	3
RS-fMRI (PA)*					
SEQUENCE NAME	GRE-EPI	TE	30	#VOLUMES	100
TYPE OF SEQUENCE	Single-shot GRE-EPI	TR	3000	FREQ. DIR	R/L
SLICE ORDER	Interleaved bottom/up	FLIP ANGLE	70	PIXEL BW	7812.5
				SLICE THICKNESS	3

*Change phase encoding direction in user CVs

To do so

Save RX > Research > Download > Display CVs

Write pepolar and set to 1